

Acrylamide					500ml		A箱
MMA			大立高分子		500ml		A箱
MMA					1L		A箱
Cinnamyl alcohol	$C_6H_5CH=CHCH_2OH$	104-54-1	景明	acros	100g		A箱
Methacrylic Acid			Janssen		500ml		A箱
DEAEA					200ml		A箱
MMA					500ml		A箱
Zink	Zn				500g		A箱
pyrrol							A箱
Norbornen	C_7H_{10}		MERCK		250g		A箱
methyl 2-bromo propionate	$(BrCH_2)_2CHCO_2CH_3$	22262-60-8	adlrich		25g		A箱
BPAF-DMAI					20g		A箱
phenylphonicdichloride		824-72-6			250ml		A箱
BPAF-DMAI					20g		A箱
Tetramethyl Tin	$C_4H_{12}Sn$	594-27-4	JANSSEN		50 gr		A箱
磷脂質		865-49-6			100ml		A箱
ACMO			KOHJIN		450g		A箱
6-Bromohexanoyl chloride	$C_6H_{10}BrClO$	22809-37-6		Aldrich	10g		A箱
BPAF-DMAI					PE瓶		A箱
BPS-DMAI					PE瓶		A箱
2-pyridine carboxyaldehyde		1121-60-4			25g		A箱

3-hydroxyphenylace		10401-11-3	lancaster		1g		A箱
BPA-DMAI					PE瓶		A箱
5-norbornene-2-endo-3-endo dimethanol			aldrich		5g		A箱
phenylmethyldichlorosilane	C6H5SiCl2CH3		tc		25g		A箱
tetramethylzinn	Sn(CH3)4				1g		A箱
Bis(tricyclohexylphosphine)benzylidene-ruthenium(IV) dichloride					1g		A箱
MAI					100ml		B箱
Iodimethan	CH3I		Merck		50mL		B箱
AIBN			500g				B箱
BPA-MAI					50ml		B箱
Methyl chloride	CH3Cl		KANIO		500gr		B箱
bis(triphenylphosphine)palladium dichloride					1g		B箱
Triisobutylaluminum	(C4H9)3Al		TCI		100mL		B箱
PTMEG-2000					1kg		B箱
CADMIUM Iodide	CdI2	7790-80-9	ACROS		250gr		B箱
DMAEA					500mL		B箱
DMAEA					500mL		B箱
dimethylacrylamide DMAA					500mL		B箱
2,2-Bipyridine	C10H8N2		merck		25g		B箱
Cholesterol	C ₂₇ H ₄₆ O		SIGMA		5g		B箱
iodine bisublime					100g		B箱

cyclopentadiene							B箱
Bromine			MERCK		250mL		C箱
Methacrylate MA			大立高分子		500g		C箱
acrylic acid					500mL		C箱
Acetyl chloride	C ₂ H ₃ ClO		acros		250mL		C箱
ethylaluminumsesquichloride	(C ₂ H ₅) ₂ AlCl	Cl ₂ AlC ₂ H ₅	aldrich		100g		C箱
methacrylic acid		79-41-4	acros		250mL		C箱
Norbornen	C ₇ H ₁₀		MERCK		250g		C箱
Norbornen	C ₇ H ₁₀		MERCK		250g		C箱
Norbornen	C ₇ H ₁₀		MERCK		250g		C箱
Copper(I) bromide	Br ₂ Cu ₂	Acros	250g				C箱
dichlorodiphenylsilane	(C ₆ H ₅) ₂ SiCl ₂	1980/10/4	aldrich		100g		C箱
Polycarbonate					500g		C箱
2-Nitrobenzoylchloride			acros		25g		C箱
Tin(II)chloride dihydrate					500g		C箱
Methyl Methacrylate MMA			大立高分子		500g		C箱
allyl alcohol					1L		C箱

3-Mercaptopropionic acid	HSCH ₂ CH ₂ CO ₂ H	107-96-0		Acros	100ml		D箱
Methacrylate MA	C ₄ H ₆ O ₂		merck		1L		D箱
1,1,1,3,3,3-hexafluoroacetone sesqui-hydrate	C ₃ F ₆ O · 1.5H ₂ O	Acros	25g				D箱
4-Nitrophenyl chloroformate,97%		7693-46-1	Aldrich				D箱
Pentaerythritol triacrylate	C=CHCO ₂ CH ₂) ₃ CCH ₂	3524-68-3	六和		5g		D箱
dichloro-n-octadecylmethysilane	CH ₃ (ch ₂) ₁₇ Si(CH ₃)Cl ₂		TCI		25ml		D箱
Benzonitrile	C ₇ H ₅ N		619	1L			側邊
styrene			merck		1L		側邊
4,4'-Diphenylmethane Diisocyanate	CH ₂ (C ₆ H ₄ NCO) ₂		merck		500g		側邊
Methacrylate MA	C ₄ H ₆ O ₂		merck		1L		側邊
Ethyl-vinylether	C ₄ H ₈ O	109-92-2	Acros		1L		側邊
Methacrylate MA			大立高分子		1L		側邊
Methyl Methacrylate MMA					1Lx2		側邊
2-EHA					1L		側邊
2-allyl phenol			acros		500g		側邊
							側邊
acrylic acid					500l		側邊
Allyl bromide	C ₃ H ₅ Br	106-95-6	Lancaster		1L		側邊
DMAPS					1kg		側邊
SPE					1kg		側邊
SPI					1kg		側邊
HEMA					500mL		側邊
pyrrolidin	C ₄ H ₉ N		merck		100mL		側邊
norbornylene,99%	C ₇ H ₁₀	448-66-8	janssen		500g		側邊
4-Vinyl-1-Cyclohexene	C ₆ H ₉ CH=CH ₂	100-40-3	merck		500mL		側邊
IPDI					100mL		側邊

移至冰箱B

Zinn (IV)-Chloride	SnCl ₄		Merck		100mL		側邊
azobis							側邊
H12MDI					1kg	林旭屏	側邊
acrylogl chloride					5g		側邊
styrene					500mL		E箱
norbornylene,99%	C ₇ H ₁₀	448-66-8	janssen		500g		E箱
2-Hydroxyethyl acrylate	C ₅ H ₈ O ₃	818-61-1		acros	1L		E箱
Allyl glycidyl ether	C ₆ H ₁₀ O ₂	106-92-3	ACROS		500g		E箱
Allyl bromide	C ₃ H ₅ Br	106-95-6	ACROS		500mL		E箱
Methyl Methacrylate MMA					500mL		E箱
Pyridin-2,6-dicarbonacide Dipicolinacide					500g		E箱
1,5-cyclooctadiene	C ₈ H ₁₂	Acros	500ml				E箱
1-hexene					250mL		E箱